

Bouffant Head Covers - (Imported)

Technical Data Sheet

KEY FEATURES

- Elastic band provides comfort and keeps cap in place on the forehead
- Comfortable lightweight SMS Nonwoven
- Spunbound/Meltblown/Spunbound

Typical Value

50g

Size

One Size Fits Most - Will comfortably accommodate 28"

Color

Blue



ORDERING INFORMATION

SKU #	DELIVERY	COVERS/PACK	PACKS/CASE	COVERS/CASE
J2BHCCN1B	Quoted Upon Request	100	10	1,000

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2020. All rights reserved.