

Technical Data Sheet

MANUFACTURER

- Intco (United States: Basic)

MATERIAL & DESIGN

- Vinyl Synthetic, latex free, powder free
- Protein free
- Color - Blue

THICKNESS

- Fingers – 0.09mm (3.54 mil)
- Palms - >= .08

LENGTH

- >= 230 mm

SIZE

- Sizes - Small, Medium, Large & Extra Large

APPROVALS AND STANDARDS

- FDA 510K Approval = K190095
- Meets ASTM Standard 5250-06 for Vinyl Gloves
- Meets ASTM Standards ASTM D6124-06 - Residual Powder
- Meets ASTM Standards: ASTM D5151-06 - Pinhole/Water Leak Test

TENSILE STRENGTH

- before/after aging 11MPA

ELONGATION

- before/after aging 300%



ORDERING INFORMATION

| | SKU # | DELIVERY | GLOVES/BOX | BOX/CASE | GLOVES/CASE |
|----|----------------|---------------------|------------|----------|-------------|
| #1 | J2VGBLSMB (SM) | | | | |
| #2 | J2VGBLMDB (MD) | Quoted Upon Request | 100 | 10 | 1,000 |
| #3 | J2VGBLLGB (LG) | | | | |
| #4 | J2VGBLXLB(XL) | | | | |

For More Information:

855.615.8633

Orders@J2MedicalSupply.com