

Disposable Nitrile Examination Glove Specs

<ul style="list-style-type: none"> • Chemo Rated • FDA 510(K): K170136 • ASTM D6319, D6978 • Powder Free • Non Sterile • Ambidextrous 	<p>Grade: Medical</p> <p>Material: NBR</p> <p>Cuff: Beaded</p> <p>Color: Blue</p> <p>Surface: Finger textured</p>	<p>Packing: 100 pcs/box,10 boxes/ctn</p> <p>Shelf Life: 5 Years</p> <p>Glove Length (mm): ≥ 230</p> <p>Weight (g): M=3.7 +/- 0.3</p> <p>Powder Level: < 2mg/glove</p>
---	--	---



Physical Dimensions

Size	Weight	Length (mm)	Palm Width (mm)	Min. Palm Thickness(mm)	Min. Finger Thickness(mm)
S	3.2 +/- 0.3	≥ 230	80 ± 10	0.06	0.08
M	3.7 +/- 0.3	≥ 230	95 ± 10	0.06	0.08
L	4.0 +/- 0.3	≥ 230	110 ± 10	0.06	0.08
XL	4.3 +/- 0.3	≥ 230	120 ± 10	0.06	0.08

Physical Properties

Force at Break (Mpa)		Elongation		Water Leak AQL	Appearance AQL
Before Aging	After Aging	Before Aging	After Aging		
≥ 14	≥ 14	≥ 500	≥ 400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥ 14	≥ 14	≥ 500	≥ 400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥ 14	≥ 14	≥ 500	≥ 400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥ 14	≥ 14	≥ 500	≥ 400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0



For More Information:

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

Chemotherapy Drug Permeation Time

	Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
1.	Acrylamide, 40%	133.3
2.	Arsenic Trioxide, 1 mg/ml (1,000 ppm)	>240
3.	Azacitidine (Vidaza), 25 mg/ml (25,000 ppm)	>240
4.	Bendamustine, 5 mg/ml (5,000 ppm)	>240
5.	Benzalkonium Chloride, 50%	>480
6.	Bleomycin, Sulfate, 15 mg/ml (15,000 ppm)	>240
7.	Blusulfan, 6 mg/ml (6,000 ppm)	>240
8.	Bortezomib (Velcade) 1 mg/ml (1,000 ppm)	>240
9.	Carboplatin, 10 mg/ml (10,000 ppm)	>240
10.	Carfilzomib, 2 mg/ml (2,000 ppm)	>240
11.	Carmustine (BCNU) (3.3 mg/ml)	9.1
12.	Cetuximab, 2 mg/ml (2,000 ppm)	240
13.	Chlorhexidine Gluconate, 4%	>480
14.	Chloroquine, 50 mg/ml (50,000 ppm)	>240
15.	Cidex OPA	66.7
16.	Cisplatin, (1.0 mg/ml)	>240
17.	Cyclophosphamide (20 mg/ml)	>240
18.	Cyclosporin A, 100 mg/ml (100,000 ppm)	>240
19.	Cytarabine, 100 mg/ml (100,000 ppm)	>240
20.	Cytovene, 10 mg/ml (10,000 ppm)	>240
21.	Dacarbazine (DTIC) (10.0 mg/ml)	>240
22.	Daunorubicin, 5 mg/ml (5,000 ppm)	>240
23.	Decitabine, 5 mg/ml (5,000 ppm)	>240
24.	Docetaxel, 10 mg/ml (10,000 ppm)	>240
25.	Doxorubicin HCl (2.0 mg/ml)	>240

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

Chemotherapy Drug Permeation Time

	Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
26.	Epirubicin HCl (Ellence), 2 mg/ml 2mg/ml (2,000 ppm)	>240
27.	Etoposide (Toposar) (20 mg/ml)	>240
28.	Fentanyl Citrate Injection, 100 mcg/2 ml	>240
29.	Fluorouracil (50 mg/ml)	>240
30.	Fluradarabine, 25 mg/ml (25,000 ppm)	>240
31.	Fulvestrant, 50 mg/ml (50,000 ppm)	>240
32.	Gemcitabine, 38 mg/ml (38,000 ppm)	>240
33.	Hydrogen Peroxide, 30%	0
34.	Idarubicin, 1 mg/ml (1,000 ppm)	>240
35.	Ifosfamide, 50 mg/ml (50,000 ppm)	>240
36.	Irinotecan, 20 mg/ml (20,000 ppm)	>240
37.	Mechlorethamine HCl, 1 mg/ml (1,000 ppm)	>240
38.	Melphalan, 5 mg/ml (5,000 ppm)	>240
39.	Mesna, 50 mg/ml (50,000 ppm)	>240
40.	Methotrexate (25 mg/ml)	>240
41.	Mitomycin C (0.5 mg/ml)	>240
42.	Mitoxantrone, 2 mg/ml (2,000 ppm)	>240
43.	Oxaliplatin, 2 mg/ml (2,000 ppm)	>240
44.	Paclitaxel (Taxol) (6.0 mg/ml)	>240
45.	Paraplatin (Carboplatin), 10 mg/ml (10,000 ppm)	>240
46.	Pemetrexed, 25 mg/ml (25,000 ppm)	>240
47.	Raltitrexed, 0.5 mg/ml (500 ppm)	>240
48.	Retrovir, 10 mg/ml (10,000 ppm)	>240
49.	Rituximab, 10 mg/ml (10,000 ppm)	>240

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

Chemotherapy Drug Permeation Time

	Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
50.	Silver Nitrate, 0.5%	>480
51.	Sodium Hypochlorite, 13%	>480
52.	Sodium Hydroxide, 40%	>480
53.	Temsirolimus, 25 mg/ml (25,000 ppm)	>240
54.	Thiotepa (10 mg/ml)	36.3
55.	Topotecan, 1 mg/ml (1,000 ppm)	>240
56.	Triclosan, 2 mg/ml (2,000 ppm)	>240
57.	Trisenox, 1 mg/ml (1,000 ppm)	>240
58.	Vinblastine, 1 mg/ml (1,000 ppm)	>240
59.	Vincristine sulfate (1 mg/ml)	>240
60.	Vinorelbine, 10 mg/ml (10,000 ppm)	>240
61.	Zoledronic Acid, 0.8 mg/ml (800 ppm)	>240

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.