

## Disposable Nitrile Examination Glove Specs

- FDA 510(K): K192333
- ASTM D6319
- Powder free
- Latex free
- Non Sterile
- Ambidextrous

**Grade:** Powder free

**Material:** NBR

**Cuff:** Beaded

**Color:** Blue

**Surface:** Finger textured

**Packing:** 100 pcs/box,10 boxes/ctn

**Shelf Life:** 5 Years

**Glove Length (mm):**  $\geq 230$

**Weight (g):** M=3.7 +/- 0.3

**Powder Level:** < 2mg/glove

## Physical Dimensions



Size	Weight	Length (mm)	Palm Width (mm)	Min. Palm Thickness(mm)	Min. Finger Thickness(mm)
S	3.4 +/- 0.3	$\geq 230$	$80 \pm 10$	0.07 +/- .02	0.09 +/- .02
M	3.9 +/- 0.3	$\geq 230$	$95 \pm 10$	0.07 +/- .02	0.09 +/- .02
L	4.2 +/- 0.3	$\geq 230$	$110 \pm 10$	0.07 +/- .02	0.09 +/- .02
XL	4.5 +/- 0.3	$\geq 230$	$120 \pm 10$	0.07 +/- .02	0.09 +/- .02

## Physical Properties

Force at Break (Mpa)		Elongation		Water Leak AQL	Appearance AQL
Before Aging	After Aging	Before Aging	After Aging		
$\geq 18$	$\geq 16$	$\geq 500$	$\geq 400$	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
$\geq 18$	$\geq 16$	$\geq 500$	$\geq 400$	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
$\geq 18$	$\geq 16$	$\geq 500$	$\geq 400$	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
$\geq 18$	$\geq 16$	$\geq 500$	$\geq 400$	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0



## For More Information:

855.615.8633

Orders@J2MedicalSupply.com

J2 Medical Supply 2021. All rights reserved.