

iHealth COVID-19 Antigen Rapid Test Technical Data Sheet

The iHealth® COVID-19 Antigen Rapid Test is a lateral flow assay intended for the qualitative detection of nucleocapsic protein antigen from SARS-CoV-2.

FEATURES

- FDA authorized 15 mins self-test
- Easy to use zero discomfort
- For ages 2 and above
- Detect current and new COVID variants
- Group testing managing via mobile app



This test is authorized for non-prescription home use with self-collected anterior nasal (nares) swab samples from individuals aged 15 years or older with symptoms of COVID-19 within the first 7 days of symptom onset.

Authorized for non-prescription home use with adult-collected nasal swab samples from individuals aged 2 years or older with symptoms of COVID-19 within the first 7 days of symptom onset.

This test is also authorized for non-prescription home use with self-collected anterior nasal (nares) swab samples from individuals aged 15 years or older, or adult collected anterior nasal swab samples from individuals aged 2 years or older, with or without symptoms or other epidemiological reasons to suspect COVID-19 when tested twice over three days with at least 24 hours (and no more than 48 hours) between tests.



Swab



Dip



Drip



Wait 15 minutes

For More Information:

855.615.8633

Orders@J2MedicalSupply.com