

Underarm Crutch

Technical Data Sheet

FEATURES

- Tube: Aluminum, Matte Color
- Handle: TPR
- Sizes - Small, Medium & Large
- Underarm Pad: TPR
- Function: Adjustable Height
- Tip: Rubber
- Height: SM (94cm-114cm), MD (114cm-134cm), LG (134cm-154cm)
- Weight Limit: 300lbs



SIZE SPECIFICATIONS:

SKU #1: J2FY925LSMC - SMALL/YOUTH

Approximate User Height: 4'6" to 5'2"

SKU #2: J2FY925LMDC - MEDIUM

Approximate User Height: 5'2" to 5'10"

SKU #3: J2FY925LLGC - LARGE

Approximate User Height: 5'10" to 6'6"

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.