

# Honeywell Disposable Respirator , DF300 Commercial Grade N95

## Technical Data Sheet

### KEY FEATURES

Engineered to promote easier breathing, the NIOSH approved DF300 Series Flatfold Disposable Respirators are made with a multilayered absorption media that repels moisture and features a smooth inner lining for added comfort.

- Designed to fit most faces comfortably and accommodate the user's facial movements
- Multilayered absorption media with a smooth inner lining for added comfort
- Latex free, comfort stretch head straps are welded for a secure fit
- Concealed, adjustable nose clip and soft foam nose cushion help provide a custom secure seal
- Humidity and moisture resistance filter media that repels moisture for a longer lasting respirator

### TECHNICAL SPECIFICATIONS

- Product weight
  - DF300 = 0.2 OZ
- Latex free components
- Comfort stretch, cloth strap
- Nose Foam: Open Cell Foam

### INTENDED USAGE

- Solid particulates & non-petroleum based liquid aerosols
- DO NOT Use for gases, vapors, oil-based mists, asbestos, paint spray, or sandblasting
- DO NOT Use for particulate concentrations exceeding 10 times PEL

### APPROVALS & STANDARDS

- NIOSH approved N95 particulate respirator (NIOSH 42 CFR 84)
  - TC-84A-8480
- Assigned Protection Factor (APF 10) per OSHA
- Meets OSHA Respiratory Protection Standard 29 CFR 1910.134

## Honeywell



### ORDERING INFORMATION

SKU #	Class Feature	Valve	Nose Clip	Strap Style	Strap Material	Box Qty	Case Qty	Case Dimensions
J2FD300HB	N95	NO	YES	Welded	Cloth	20	18	20.5"L x15.8"W x12"H

### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2020. All rights reserved.