

Protective Coverall

Technical Data Sheet

PROTECTION

The Coverall provides protection against dry particles and liquids in multi-use environments with an AAMI Level 4 rating.

PROTECTION

The Coverall is made of advanced PP+PE material that offers a combination of comfort and protection for demanding situations. This durable protective coverall meets AAMI PB70 Level 4 testing requirements.

KEY FEATURES

- Designed for freedom of movement with extra material where you need it most.
- Long sleeves with snug fit elastic wrists.
- Full length zipper in front and elastic gatherings in back for a superior fit.
- Strong, ultrasonic sealed seams are designed to be durable and offer complete protection.



PRODUCT SPECIFICATIONS

Product Codes	CVRL4-HDNB-3R, CVRL4-
Material	Heavy Weight PP+PE
Fabric Weight	70G
Packaging	Individually Packaged, 25/case
Origin	Made in Cambodia, TAA Compliant

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.