

## Disposable Lab Coat

### Technical Data Sheet

---

Lab Coat can be used in a variety of healthcare settings and provides protection against dry particles and low levels of fluid.

#### PROTECTION

Made of advanced, light weight spunbond polypropylene material and offers protections against fluids and spray, while maintaining breathability and improved air circulation.

#### GOWN FEATURES

Generously sized, knee-length lab coat comes rounded knit collar and knit cuffs. Five button front closure affords a comfortable a strong and durable fit and offers complete protection. Low-linting material ensures a clean environment. Does not contain components made from rubber latex.



---

#### PRODUCT SPECIFICATIONS

<b>Material</b>	Spunbond Polypropylene
<b>Fabric Weight</b>	35 gsm
<b>Packaging</b>	10/bag, 100/case, 24 cases/pallet
<b>Case Dimensions</b>	20.5x12.6x21.3in, GW 22.5 lbs

#### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.