

J2-100 Advanced

Disposable Nitrile Examination Glove Specs

<ul style="list-style-type: none"> • ASTM D6978 • Powder Free • Non Sterile • Ambidextrous • Chemo Rated 	Grade: Medical Grade	Packing: 100 pcs/box,10 boxes/ctn
	Material: NBR	Shelf Life: 5 Years
	Cuff: Rolled Beaded	Glove Length (mm): ≥230
	Color: Blue	Weight (g): M=3.7 +/- 0.3
	Surface: Finger textured	Powder Level: < 2mg/glove

Physical Dimensions



Size	Weight	Length (mm)	Palm Width (mm)	Min. Palm Thickness(mm)	Min. Finger Thickness(mm)
S	4.0 +/- 0.3	≥240	80 ± 10	0.08 +/- 0.02	0.11+/0.01
M	4.5 +/- 0.3	≥240	95 ± 10	0.08 +/- 0.02	0.11+/0.01
L	5.0 +/- 0.3	≥240	110 ± 10	0.08 +/- 0.02	0.11+/0.01
XL	5.5 +/- 0.3	≥240	120 ± 10	0.08 +/- 0.02	0.11+/0.01

Physical Properties

Force at Break (Mpa)		Elongation		Water Leak AQL	Appearance AQL
Before Aging	After Aging	Before Aging	After Aging		
≥15	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥15	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥15	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥15	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0



For More Information:

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.