

# Face Shields

## Technical Data Sheet

### KEY FEATURES

- Lightweight, anti-fog protection from splashing or spraying
- Foam Headband with Elastic Strap
- Can be used in conjunction with a N95 mask, goggles, glasses or be worn alone without fogging up
- Made from polyethylene terephthalate (PET)

<b>Color</b>	Clear Shield / Blue Headband
<b>Dimensions</b>	9 Inch Length Shield / 1 X 1 X 9 Inch Foam Strip
<b>Fastening Type</b>	Foam Headband / Elastic Strap
<b>Lens Tint</b>	Clear Tint
<b>Size</b>	One Size Fits Most
<b>Sterility</b>	NonSterile
<b>Style</b>	Full Length
<b>Style</b>	Anti-fog



### APPROVALS AND STANDARDS

- Compliance to ANSI Z87.1 -2015 Occupational and Educational Personal Eye and Face Protection Devices
- Clause 5.1.2 - Luminous Transmission
- Clause 8.1.2 - Droplet and Splash Hazard

### ORDERING INFORMATION

SKU #	DELIVERY	SHIELDS/PACK	PACKS/CASE	SHIELDS/CASE
J2FSMG01C	Quoted Upon Request	12	20	240

#### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.