


## Moldex 1500 N95 Series Technical Data Sheet



1510

1511

1512

1513

1517

### ITEM GENERAL DESCRIPTION:

N95 Healthcare Particulate Respirator and Surgical Mask

MODEL NUMBER:

1510, 1511, 1512, 1513, 1517

THIS REVISION ISSUE DATE:

2/05/15

REVISION LETTER: 2

CDR: 15-015

### MATERIALS USED:

1. EVA/Polyethylene
2. Polypropylene
3. Polyethylene
4. Tin Plated Steel Wire
5. TPE

### PARTS WERE USED:

1. Dura-Mesh® Shell
2. Filter Media
3. Softspun® Inner Liner
4. Nose Flange
5. Staples
6. Straps

### SPECIFIC PRODUCT FUNCTION

Respirator is intended as a particulate filter at protection level N95 and surgical mask with limitations as noted on the container. NIOSH certified for 1510 per approval number TC-84A-5171 and NIOSH certified for 1511, 1512, 1513 and 1517 per approval number TC-84A-0013.

### FIRE DATA

Extinguishing Media – Water Spray, CO2, Foam, Dry chemical

### ENVIRONMENTAL INFORMATION

Recycle if possible. Product is landfill disposable but when used in healthcare environments this product may have to be treated as a biohazard. For all disposal, consult with state and local authorities.

### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.