


### Technical Data Sheet

#### KEY FEATURES

- SMS Nonwoven - 35gsm
- Over the head with thumb loop
- Perforated back for easy removal
- Product Color: Yellow
- Available in two sizes - SM/MD and LG/XL
- Superior comfort while providing protection from light to moderate fluid exposure
- Product must be ordered in full case increments, as we do not break cases on direct sales



#### AAMI PB70 Guidelines

	LEVEL 1	LEVEL 2		LEVEL 3	
Test	AATCC 42:2000	AATCC 42:2000	AATCC 127:1998	AATCC 42:2000	AATCC 127:1998
Requirements at 4% AQL	Water Impact ≤ 4.5 g	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 20 cm	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 50 cm
Anticipated Risk of Exposure	Low	Moderate		Moderate to High	


For additional information on guidelines on gowns, please refer to the following FDA link - <https://www.fda.gov/medical-devices/personal-protective-equipment-infection-control/medical-gowns>

#### For More Information:


855.615.8633

Orders@J2MedicalSupply.com

## DESIGN: Sunshine SM-MD



## DESIGN: Sunshine LG-XL



### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.

## INTENDED USAGE

- Single-use gown to protect the wearer from the transfer of microorganisms and body fluids in low or minimal risk patient isolation situations.
- For use in patient care areas and other points of interaction where there is a need for barrier protection.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.



## ORDERING INFORMATION

	SKU #	DELIVERY	BOX/CASE COUNT
#1	J2IG2SSMC (SM/MD)	Quoted Upon Request	100 units per case
#2	J2IG2SLXC (LG/XL)		

### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.