

AAMI Level 4 Sterile Surgical Gown

Technical Data Sheet

KEY FEATURES

- Knit cuffs and side tie
- Full back coverage with hook and loop neck closure
- Heat-sealed, ultrasonic and sewn seams
- Low linting construction
- Does not contain components made from rubber latex
- Product Color: Blue

PROTECTION

This sterile level 4 reinforced surgical gown provides great protection in a lightweight material with excellent moisture evaporation. Constructed with reinforced chest and sleeves. Protects against microorganisms, liquids and particles in operating room settings.



RATING LEVELS

	Good	Better	Best
Protection			●
Comfort			●

PRODUCT SPECIFICATIONS

Material	Medium weight SMMS
Fabric Weight	50 gsm
Packed	Ships with disposable surgical towel
510k No.	K121998
Origin	Made in China

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.