

## Underarm Crutch Technical Data Sheet

### FEATURES

- Tube: Aluminum, Matte Color
- Handle: TPR
- Sizes - Small, Medium & Large
- Underarm Pad: TPR
- Function: Adjustable Height
- Tip: Rubber
- Height: SM (94cm-114cm), MD (114cm-134cm), LG (134cm-154cm)
- Weight Limit: 300lbs
- Country of Origin: China



### SIZE SPECIFICATIONS:

**SKU #1:** [J2FY925LSMC - SMALL/YOUTH](#)

Approximate User Height: 4'6" to 5'2"

**SKU #2:** [J2FY925LMDC - MEDIUM](#)

Approximate User Height: 5'2" to 5'10"

**SKU #3:** [J2FY925LLGC - LARGE](#)

Approximate User Height: 5'10" to 6'6"

### For More Information:

855.615.8633

[Orders@J2MedicalSupply.com](mailto:Orders@J2MedicalSupply.com)

© J2 Medical Supply 2021. All rights reserved.