

Clear Vinyl Examination Gloves

Technical Data Sheet

COUNTRY OF ORIGIN

- China

MATERIAL & DESIGN

- Vinyl Synthetic, latex free, powder free
- Protein free
- Color - Clear

THICKNESS

- 4 ml
- Fingers – 0.09mm (3.54 mil)
- Palms - 0.07mm (2.76 mil)

LENGTH

- ≥ 240 mm

SIZE

- Sizes - Small, Medium, Large & Extra Large

APPROVALS AND STANDARDS

- FDA 510K Approval = K190095
- Meets ASTM Standard 5250-06 for Vinyl Gloves
- Meets ASTM Standards ASTM D6124-06 - Residual Powder
- Meets ASTM Standards: ASTM D5151-06 - Pinhole/Water Leak Test

TENSILE STRENGTH

- Before/after aging 11MPA

ELONGATION

- Before/after aging 300%



ORDERING INFORMATION

	SKU #	DELIVERY	GLOVES/BOX	BOX/CASE	GLOVES/CASE
#1	J2VG4MSMC (SM)	Quoted Upon Request	100	10	1,000
#2	J2VG4MMDC (MD)				
#3	J2VG4MLGC (LG)				
#4	J2VG4MXLC (XL)				

For More Information:

855.615.8633

Orders@J2MedicalSupply.com