

Protective Coverall

Technical Data Sheet

PROTECTION

The Coverall provides protection against dry particles and liquids in multi-use environments with an AAMI Level 3 rating.

The Coverall is made of advanced PP+PE material that offers a combination of comfort and protection for demanding situations. This durable protective coverall meets AAMI PB70 Level 3 testing requirements.

KEY FEATURES

- Designed for freedom of movement with extra material where you need it most.
- Long sleeves with snug fit elastic wrists.
- Full length zipper in front and elastic gatherings in back for a superior fit.
- Strong, ultrasonic sealed seams are designed to be durable and offer complete protection.



ORDERING INFORMATION

	SKU #	GTIN	DELIVERY	BOX/CASE COUNT
#1	J2L4COVRC (LG)	00850020592883	Quoted Upon Request	25 units per case
#2	J2L4COVXC (XL)	00850020592890		
#3	J2L4COV2C (2XL)	00850020592906		
#4	J2L4COV3C (3XL)	00850020592913		

PRODUCT SPECIFICATIONS

Product Codes	CVRL4-HDNC-3R, CVRL4-
Material	Heavy Weight PP+PE
Fabric Weight	70G
Packaging	Individually Packaged, 25/case
Origin	Made in Cambodia, TAA Compliant

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.