

Honeywell Disposable Respirator, DC365 Surgical N95

Technical Data Sheet

KEY FEATURES

The Honeywell DC365 Surgical N95 Respirator is an easy-to-don and doff, high fluid-resistant surgical respirator designed for clinician satisfaction across the spectrum of care.

- Fluid-resistant at 160mmHg• 160mmHg
- Designed with comfort in mind• mind
- Dye-free and non-sensitizing• sensitizing
- Materials feel smooth and soft next to the skin• skin
- Flexible, light-weight construction conforms gently to the face, minimizing pressure pointspressure points
- Cup-shape with spacious breathing chamber
- Durable, dual-welded straps are easy to don and doff and provide a secure and comfortable fit
- Adjustable nosepiece for a customizable fit
- Multi-layer absorption media designed for lower breathing resistance• resistance
- One size

TECHNICAL SPECIFICATIONS

- NIOSH-certified Surgical N95 (NIOSH 42 CFR 84)
- FDA registered
- Fluid-resistant ASTM F1862 160 mmHg

INTENDED USAGE

Honeywell delivers high-performance solutions that are innovative and designed for exceptional fit. We've spent decades researching and testing new materials to develop the best PPE for the broadest range of user needs. Our experience, coupled with our insatiable desire to adapt and help protect workers in every environment, means healthcare professionals can rely on Honeywell's history of quality and innovation to help keep them protected while on the job.

Honeywell



ORDERING INFORMATION

SKU #	Class Feature	Box Qty	Case Qty	Units Case
J2DC365HC	N95	20	10	200

For More Information:

855.615.8633

Orders@J2MedicalSupply.com