

## J2-100 Premium Nitrile

### DISPOSABLE NITRILE EXAMINATION GLOVE SPECS

- FDA 510(K): K170136
- ASTM D6978
- Powder Free
- Non Sterile
- Ambidextrous
- Chemo Rated

**Grade:** Medical

**Material:** NBR

**Cuff:** Beaded

**Color:** Blue

**Surface:** Finger textured

**Packing:** 100 pcs/box,10 boxes/ctn

**Shelf Life:** 5 Years

**Glove Length (mm):** ≥230

**Weight (g):** M=3.7 +/- 0.3

**Powder Level:** < 2mg/glove



### PHYSICAL DIMENSIONS

Size	Weight	Length (mm)	Palm Width (mm)	Min. Palm Thickness(mm)	Min. Finger Thickness(mm)
S	3.2 +/- 0.3	≥230	80 ± 10	0.06	0.08
M	3.7 +/- 0.3	≥230	95 ± 10	0.06	0.08
L	4.0 +/- 0.3	≥230	110 ± 10	0.06	0.08
XL	4.3 +/- 0.3	≥230	120 ± 10	0.06	0.08

### PHYSICAL PROPERTIES

Force at Break (Mpa)		Elongation		Water Leak AQL	Appearance AQL
Before Aging	After Aging	Before Aging	After Aging		
≥14	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥14	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥14	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥14	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0



**For More Information:**

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

## Technical Data Sheet - (Continued)

Chemotherapy Drug Permeation Time

	<b>Chemotherapy Drug and Concentration</b>	<b>Minimum Breakthrough Detection Time (Minutes)</b>
1.	Acrylamide, 40%	133.3
2.	Arsenic Trioxide, 1 mg/ml (1,000 ppm)	>240
3.	Azacitidine (Vidaza), 25 mg/ml (25,000 ppm)	>240
4.	Bendamustine, 5 mg/ml (5,000 ppm)	>240
5.	Benzalkonium Chloride, 50%	>480
6.	Bleomycin, Sulfate, 15 mg/ml (15,000 ppm)	>240
7.	Blusulfan, 6 mg/ml (6,000 ppm)	>240
8.	Bortezomib (Velcade) 1 mg/ml (1,000 ppm)	>240
9.	Carboplatin, 10 mg/ml (10,000 ppm)	>240
10.	Carfilzomib, 2 mg/ml (2,000 ppm)	>240
11.	Carmustine (BCNU) (3.3 mg/ml)	9.1
12.	Cetuximab, 2 mg/ml (2,000 ppm)	240
13.	Chlorhexidine Gluconate, 4%	>480
14.	Chloroquine, 50 mg/ml (50,000 ppm)	>240
15.	Cidex OPA	66.7
16.	Cisplatin, (1.0 mg/ml)	>240
17.	Cyclophosphamide (20 mg/ml)	>240
18.	Cyclosporin A, 100 mg/ml (100,000 ppm)	>240
19.	Cytarabine, 100 mg/ml (100,000 ppm)	>240
20.	Cytovene, 10 mg/ml (10,000 ppm)	>240
21.	Dacarbazine (DTIC) (10.0 mg/ml)	>240
22.	Daunorubicin, 5 mg/ml (5,000 ppm)	>240
23.	Decitabine, 5 mg/ml (5,000 ppm)	>240
24.	Docetaxel, 10 mg/ml (10,000 ppm)	>240
25.	Doxorubicin HCl (2.0 mg/ml)	>240

**For More Information:**

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

## Technical Data Sheet - (Continued)

Chemotherapy Drug Permeation Time

	<b>Chemotherapy Drug and Concentration</b>	<b>Minimum Breakthrough Detection Time (Minutes)</b>
26.	Epirubicin HCl (Ellence), 2 mg/ml 2mg/ml (2,000 ppm)	>240
27.	Etoposide (Toposar) (20 mg/ml)	>240
28.	Fentanyl Citrate Injection, 100 mcg/2 ml	>240
29.	Fluorouracil (50 mg/ml)	>240
30.	Fluradarabine, 25 mg/ml (25,000 ppm)	>240
31.	Fulvestrant, 50 mg/ml (50,000 ppm)	>240
32.	Gemcitabine, 38 mg/ml (38,000 ppm)	>240
33.	Hydrogen Peroxide, 30%	0
34.	Idarubicin, 1 mg/ml (1,000 ppm)	>240
35.	Ifosfamide, 50 mg/ml (50,000 ppm)	>240
36.	Irinotecan, 20 mg/ml (20,000 ppm)	>240
37.	Mechlorethamine HCl, 1 mg/ml (1,000 ppm)	>240
38.	Melphalan, 5 mg/ml (5,000 ppm)	>240
39.	Mesna, 50 mg/ml (50,000 ppm)	>240
40.	Methotrexate (25 mg/ml)	>240
41.	Mitomycin C (0.5 mg/ml)	>240
42.	Mitoxantrone, 2 mg/ml (2,000 ppm)	>240
43.	Oxaliplatin, 2 mg/ml (2,000 ppm)	>240
44.	Paclitaxel (Taxol) (6.0 mg/ml)	>240
45.	Paraplatin (Carboplatin), 10 mg/ml (10,000 ppm)	>240
46.	Pemetrexed, 25 mg/ml (25,000 ppm)	>240
47.	Raltitrexed, 0.5 mg/ml (500 ppm)	>240
48.	Retrovir, 10 mg/ml (10,000 ppm)	>240
49.	Rituximab, 10 mg/ml (10,000 ppm)	>240

**For More Information:**

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

## Technical Data Sheet - (Continued)

Chemotherapy Drug Permeation Time

	<b>Chemotherapy Drug and Concentration</b>	<b>Minimum Breakthrough Detection Time (Minutes)</b>
50.	Silver Nitrate, 0.5%	>480
51.	Sodium Hypochlorite, 13%	>480
52.	Sodium Hydroxide, 40%	>480
53.	Temsirolimus, 25 mg/ml (25,000 ppm)	>240
54.	Thiotepa (10 mg/ml)	36.3
55.	Topotecan, 1 mg/ml (1,000 ppm)	>240
56.	Triclosan, 2 mg/ml (2,000 ppm)	>240
57.	Trisenox, 1 mg/ml (1,000 ppm)	>240
58.	Vinblastine, 1 mg/ml (1,000 ppm)	>240
59.	Vincristine sulfate (1 mg/ml)	>240
60.	Vinorelbine, 10 mg/ml (10,000 ppm)	>240
61.	Zoledronic Acid, 0.8 mg/ml (800 ppm)	>240

**For More Information:**

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

**COUNTRY OF ORIGIN**

- Made in China

**ORDERING INFORMATION - CASES**

	SKU #	GTIN	DELIVERY	BOX/CASE COUNT
#1	JVBBCHSMC (SM)	00850020592555	Quoted Upon Request	1000 units per case
#2	JVBBCHMDC (MD)	00850020592562		
#3	JVBBCHLGC (LG)	00850020592579		
#4	JVBBCHXLC (XL)	00850020592586		

**ORDERING INFORMATION - BOXES**

	SKU #	GTIN	DELIVERY	BOX/CASE COUNT
#1	JVBBCHSMB (SM)	00850020592593	Quoted Upon Request	Product must be ordered in full case increments (10 boxes in each)
#2	JVBBCHMDB (MD)	00850020592609		
#3	JVBBCHLGB (LG)	00850020592616		
#4	JVBBCHXLB (XL)	00850020592623		

**For More Information:**

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.