

KEY FEATURES

- NIOSH approved N95 Rating
- This N95 is FDA approved as surgical mask
- 4 layers design
- Adjustable nose clip significantly reduces air penetration
- Ultrasonically welded head strap
- Comfortable urethane foam nose-piece to reduce leakage
- Sizes - Standard and Small

MATERIAL COMPOSITION

- Outer and inner layers -Polyethylene
- Filter (two layers)-Melt-blown Nonwoven Polypropylene
- Straps -Polyester Elastic Cord
- Nose clip -Aluminum
- Nose pad -Urethane Foam
- This mask contains no components made from natural rubber latex
- Approximate weight of product: 8.5 grams / 0.3 oz



COUNTRY OF ORIGIN

- Made in China

INTENDED USAGE

- Protective mask is intended to be worn by general public or healthcare personnel and cover the user's nose and mouth.
- The mask is used for protection against particles, such as those from minerals, coal, iron, ore, flour and certain other substances as well as liquid or non-oil based particles from sprays that do not also emit oil aerosols or vapors.

DO NOT USE FOR

- The mask is not for gases, vapors and oil aerosols.
- The mask is not for particulate concentration exceeding 10X PEL/OEL.
- The respirator is not guaranteed to protect against the risk of contracting a disease or infection.

TIME USE LIMITATION

- Replace the mask when it becomes dirty, damaged, or difficult to breathe through.

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.

N95 Respirator Masks, Makrite 9500

Technical Data Sheet - (Continued)

APPROVALS AND STANDARDS

- NIOSH Approved N95 Particulate Respirator
- Meets NIOSH 42 CFR 84 N95 requirements for a minimum 95% filtration efficiency
- NIOSH Approval Number: TC-84A-5411
- FDA 510(K) Exemption

SHELF LIFE AND STORAGE

- 3 years from the date of manufacture
- Useby date on box in MM/YYYY format
- Store masks in the original packaging, away from contaminated areas, dust, sunlight,extreme temperatures, excessive moisture, and damaging chemicals
- Store in temperatures between -4°F (-20°C) and +86°F (+30°C) and not exceeding80% humidity



WARNINGS

- The high protection masks are not intended for reuse and shall be discarded after each procedure. Check your mask's production date. Masks expire two years after the production date, after which they may not provide full protection. Always follow the instructions for use to ensure mask is properly sealed. Not recommended for children under 3.

ORDERING INFORMATION

	SKU #	GTIN	DELIVERY	MASKS/BOX	BOXES/CASE	MASKS/CASE
#1	J2N95RM7C (Standard)	00850020592319	Quoted Upon Request	20	12	240
#2	J2N95RM6C (Small)	00850020592326				

Listed on CDC NIOSH approved N95 Particulate Filter Face Piece Respirator List - Approval No. 84A-5411 (Standard)

Listed on CDC NIOSH approved N95 Particulate Filter Face Piece Respirator List - Approval No. 84A-5463 (Small)

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.

J2 | MEDICAL
SUPPLY
J2MedicalSupply.com