


Technical Data Sheet

KEY FEATURES

- Full back coverage
- Elastic on the wrist for comfortable fit
- Provides comfort while providing protection from moderate to high fluid exposure
- 50gsm Polypropylene Non-woven
- Sizes - One Size Fits Most
- Product Color: Blue



AAMI PB70 Guidelines

| | LEVEL 1 | LEVEL 2 | | LEVEL 3 | |
|------------------------------|-------------------------|-------------------------|---------------------------------|-------------------------|---------------------------------|
| Test | AATCC 42:2000 | AATCC 42:2000 | AATCC 127:1998 | AATCC 42:2000 | AATCC 127:1998 |
| Requirements at 4% AQL | Water Impact ≤ 4.5 g | Spray Impact ≤ 1.0 g | Hydrostatic Pressure ≥ 20 cm | Spray Impact ≤ 1.0 g | Hydrostatic Pressure ≥ 50 cm |
| Anticipated Risk of Exposure | Low | Moderate | | Moderate to High | |

For additional information on guidelines on gowns, please refer to the following FDA link - <https://www.fda.gov/medical-devices/personal-protective-equipment-infection-control/medical-gowns>

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

INTENDED USAGE

- Single-use gown to protect the wearer from the transfer of microorganisms and body fluids in low or minimal risk patient isolation situations.
- For use in patient care areas and other points of interaction where there is a need for barrier protection.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.

COUNTRY OF ORIGIN

- Made in China



ORDERING INFORMATION

| SKU # | GTIN | DELIVERY | BOX/CASE COUNT |
|-----------|----------------|---------------------|--------------------|
| J2IG3SSMC | 00850020592357 | Quoted Upon Request | 100 units per case |

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.