

Technical Data Sheet

KEY FEATURES

- Single-layer polyethylene
- Reliable protection from transfer of microorganisms and body fluids
- Over the head with thumb loop
- Waist and neck tie closures
- Lightweight, breathable and comfortable
- Latex Free
- Sizes - One Size Fits Most
- Product Color: Blue



AAMI PB70 Guidelines

	LEVEL 1	LEVEL 2		LEVEL 3	
Test	AATCC 42:2000	AATCC 42:2000	AATCC 127:1998	AATCC 42:2000	AATCC 127:1998
Requirements at 4% AQL	Water Impact ≤ 4.5 g	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 20 cm	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 50 cm
Anticipated Risk of Exposure	Low	Moderate		Moderate to High	

For additional information on guidelines on gowns, please refer to the following FDA link - <https://www.fda.gov/medical-devices/personal-protective-equipment-infection-control/medical-gowns>

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

INTENDED USAGE

- Single-use gown to protect the wearer from the transfer of microorganisms and body fluids in low or minimal risk patient isolation situations.
- For use in patient care areas and other points of interaction where there is a need for barrier protection.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.

COUNTRY OF ORIGIN

- Made in China

GTIN

- 00850020592029



ORDERING INFORMATION

SKU #	GTIN	DELIVERY	BOX/CASE COUNT
J2IG20CNC	00850020592029	Quoted Upon Request	100 units per case

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2021. All rights reserved.