

Disposable Nitrile Examination Glove Specs

<ul style="list-style-type: none"> • ASTM D6319 • Powder free • Latex free • Non Sterile • Ambidextrous • Tested for use with Fentanyl 	<p>Grade: Medical Grade</p> <p>Material: NBR</p> <p>Cuff: Rolled-Beaded</p> <p>Color: Black</p> <p>Surface: Finger textured</p>	<p>Packing: 100 pcs/box,10 boxes/ctn</p> <p>Shelf Life: 5 Years</p> <p>Glove Length (mm): >=240</p> <p>Weight (g): M>=5.5g</p> <p>Powder Level: ≤2mg/glove</p> <p>Country of Origin: China</p>
--	--	--

Physical Dimensions



Size	Weight	Length (mm)	Palm Width (mm)	Min. Palm Thickness(mm)	Min. Finger Thickness(mm)
S	>=5.0	>=240	80±10	0.11+/0.03	0.14+/0.03
M	>=5.5	>=240	95±10	0.11+/0.03	0.14+/0.03
L	>=6.0	>=240	110±10	0.11+/0.03	0.14+/0.03
XL	>=6.5	>=240	120±10	0.11+/0.03	0.14+/0.03

Physical Properties

Force at Break (Mpa)		Elongation		Water Leak AQL	Appearance AQL
Before Aging	After Aging	Before Aging	After Aging		
>=15	>=14	>=500	>=400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
>=15	>=14	>=500	>=400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
>=15	>=14	>=500	>=400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
>=15	>=14	>=500	>=400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0



For More Information:

855.615.8633

Orders@J2MedicalSupply.com

J2 Medical Supply 2021. All rights reserved.