

Nitrile Examination Glove Technical Data Sheet J2-100 Premium Nitrile

DISPOSABLE NITRILE EXAMINATION GLOVE SPECS

- FDA 510(K): K170136
- ASTM D6978
- Powder Free
- Non Sterile
- Ambidextrous
- Chemo Rated

Grade: Medical

Material: NBR

Cuff: Beaded

Color: Blue

Surface: Finger textured

Packing: 100 pcs/box, 10 boxes/ctn

Shelf Life: 5 Years

Glove Length (mm): ≥ 230

Weight (g): M=3.7 +/- 0.3

Powder Level: < 2mg/glove



PHYSICAL DIMENSIONS

Size	Weight	Length (mm)	Palm Width (mm)	Min. Palm Thickness(mm)	Min. Finger Thickness(mm)
S	3.2 +/- 0.3	≥ 230	80 ± 10	0.06	0.08
M	3.7 +/- 0.3	≥ 230	95 ± 10	0.06	0.08
L	4.0 +/- 0.3	≥ 230	110 ± 10	0.06	0.08
XL	4.3 +/- 0.3	≥ 230	120 ± 10	0.06	0.08

PHYSICAL PROPERTIES

Force at Break (Mpa)		Elongation		Water Leak AQL	Appearance AQL
Before Aging	After Aging	Before Aging	After Aging		
≥ 14	≥ 14	≥ 500	≥ 400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥ 14	≥ 14	≥ 500	≥ 400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥ 14	≥ 14	≥ 500	≥ 400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0
≥ 14	≥ 14	≥ 500	≥ 400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0



For More Information:

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

Technical Data Sheet - (Continued)

Chemotherapy Drug Permeation Time

	Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
1.	Arsenic Trioxide (1 mg/ml)	>240
2.	Azacitidine (Vidaza) (25 mg/ml)	>240
3.	Bendamustine (5 mg/ml)	>240
4.	Bleomycin Sulfate (15 mg/ml)	>240
5.	Bortezomib (Velcade) (1 mg/ml)	>240
6.	Busulfan (6 mg/ml)	>240
7.	Carboplatin (10 mg/ml)	>240
8.	Carfilzomib (2 mg/ml)	>240
9.	Carmustine (3.3mg/ml)	34.7
10.	Cetuximab (2 mg/ml)	>240
11.	Chloroquine (50 mg/ml)	>240
12.	Cisplatin (1 mg/ml)	>240
13.	Cladribine (1 mg/ml)	>240
14.	Cyclophosphamide (20 mg/ml)	>240
15.	Cyclosporine-Organ Rejection Drug (100 mg/ml)	>240
16.	Cytarabine (100 mg/ml)	>240
17.	Cytovene (10 mg/ml)	>240
18.	Dacarbazine (10 mg/ml)	>240
19.	Daunorubicin (5 mg/ml)	>240
20.	Decitabine (5 mg/ml)	>240
21.	Docetaxel (10 mg/ml)	>240
22.	Doxorubicin HCl (2 mg/ml)	>240
23.	Epirubicin HCl (Ellence) (2 mg/ml)	>240
24.	Etoposide (20 mg/ml)	>240
25.	Fludarabine (25 mg/ml)	>240

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

Technical Data Sheet - (Continued)

Chemotherapy Drug Permeation Time

	Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
26.	Fluorouracil (50 mg/ml)	>240
27.	Fulvestrant (50 mg/ml)	>240
28.	Gemcitabine (38 mg/ml)	>240
29.	Idarubicin (1 mg/ml)	>240
30.	Ifosfamide (50 mg/ml)	>240
31.	Irinotecan (20 mg/ml)	>240
32.	Melphalan (5 mg/ml)	>240
33.	Mesna (50 mg/ml)	>240
34.	Mechlorethamine HCL (1mg/ml)	>240
35.	Methotrexate (25 mg/ml)	>240
36.	Mitomycin C (0.5 mg/ml)	>240
37.	Mitoxantrone (2 mg/ml)	>240
38.	Oxaliplatin (2 mg/ml)	>240
39.	Paclitaxel (6 mg/ml)	>240
40.	Paraplatin (10 mg/ml)	>240
41.	Pemetrexed (25 mg/ml)	>240
42.	Raltitrexed (0.5 mg/ml)	>240
43.	Retrovir (10 mg/ml)	>240
44.	Rituximab (10 mg/ml)	>240
45.	Temsirolimus (25 mg/ml)	>240
46.	ThioTEPA (10 mg/ml)	47.8
47.	Topotecan (1 mg/ml)	>240
48.	Triclosan (2 mg/ml)	>240
49.	Trisenox (1 mg/ml)	>240

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

Technical Data Sheet - (Continued)

Chemotherapy Drug Permeation Time

	Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
50.	Vinblastine (1 mg/ml)	>240
51.	Vincristine Sulfate (1 mg/ml)	>240
52.	Vinorelbine (10 mg/ml)	>240
53.	Zoledronic Acid (0.8 mg/ml)	>240

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

Technical Data Sheet - (Continued)

COUNTRY OF ORIGIN

- Made in China

ORDERING INFORMATION - CASES

	SKU #	GTIN	DELIVERY	BOX/CASE COUNT
#1	JVBBCHSMC (SM)	00850020592555	Quoted Upon Request	1000 units per case
#2	JVBBCHMDC (MD)	00850020592562		
#3	JVBBCHLGC (LG)	00850020592579		
#4	JVBBCHXLC (XL)	00850020592586		

ORDERING INFORMATION - BOXES

	SKU #	GTIN	DELIVERY	BOX/CASE COUNT
#1	JVBBCHSMB (SM)	00850020592432	Quoted Upon Request	Product must be ordered in full case increments (10 boxes in each)
#2	JVBBCHMDB (MD)	00850020592449		
#3	JVBBCHLGB (LG)	00850020592456		
#4	JVBBCHXLB (XL)	00850020592463		

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.