

Nitrile Examination Glove Technical Data Sheet J2-100 Advanced

Disposable Nitrile Examination Glove Specs

- Chemo Rated
- FDA 510(K): K170136
- ASTM D6319, D6978
- Powder Free
- Non Sterile
- Ambidextrous

Grade: Medical

Material: NBR

Cuff: Rolled-Beaded

Color: Blue

Surface: Finger textured

Water Leak AQL: 1.5

Appearance AQL: Critical AQL1.0 Major AQL2.5 Minor 4.0

Packing: 100 EA/BX | 10 BX/CA
200 EA/BX* | 10 BX/CA

Shelf Life: 5 Years

Glove Length (mm): ≥ 230

Weight (g): M=3.7 +/- 0.3

Powder Level: < 2mg/glove

*Size XL = 180 EA/BX

Physical Dimensions



Size	Weight	Length (mm)	Palm Width (mm)	Min. Palm Thickness (mm)	Min. Finger Thickness (mm)
Small	3.2+/-0.3	≥ 230	80+/-10	0.06	0.08
Medium	3.7+/-0.3	≥ 230	95+/-10	0.06	0.08
Large	4.0+/-0.3	≥ 230	110+/-10	0.06	0.08
X-Large	4.3+/-0.3	≥ 230	120+/-10	0.06	0.08

Physical Properties

Size	Force at Break (Mpa)		Elongation	
	Before Aging	After Aging	Before Aging	After Aging
Small	≥ 14	≥ 14	≥ 500	≥ 400
Medium	≥ 14	≥ 14	≥ 500	≥ 400
Large	≥ 14	≥ 14	≥ 500	≥ 400
X-Large	≥ 14	≥ 14	≥ 500	≥ 400



COUNTRY OF ORIGIN: Made in China

For More Information:

855.615.8633

Info@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.

Chemotherapy Drug Permeation Time

	Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
1	Arsenic Trioxide (1 mg/ml)	>240
2	Azacitidine (Vidaza) (25 mg/ml)	>240
3	Bendamustine (5 mg/ml)	>240
4	Bleomycin Sulfate (15 mg/ml)	>240
5	Bortezomib (Velcade) (1 mg/ml)	>240
6	Busulfan (6 mg/ml)	>240
7	Carboplatin (10 mg/ml)	>240
8	Carfilzomib (2 mg/ml)	>240
9	Carmustine (3.3mg/ml)	34.7
10	Cetuximab (2 mg/ml)	>240
11	Chloroquine (50 mg/ml)	>240
12	Cisplatin (1 mg/ml)	>240
13	Cladribine (1 mg/ml)	>240
14	Cyclophosphamide (20 mg/ml)	>240
15	Cytarabine (100 mg/ml)	>240
16	Cytovene (10 mg/ml)	>240
17	Dacarbazine (10 mg/ml)	>240
18	Daunorubicin (5 mg/ml)	>240
19	Decitabine (5 mg/ml)	>240
20	Docetaxel (10 mg/ml)	>240
21	Doxorubicin HCl (2 mg/ml)	>240
22	Epirubicin HCl (Ellence) (2 mg/ml)	>240
23	Etoposide (20 mg/ml)	>240
24	Fludarabine (25 mg/ml)	>240
25	Fluorouracil (50 mg/ml)	>240
26	Fulvestrant (50 mg/ml)	>240
27	Gemcitabine (38 mg/ml)	>240
28	Idarubicin (1 mg/ml)	>240
29	Ifosfamide (50 mg/ml)	>240
30	Irinotecan (20 mg/ml)	>240
31	Melphalan (5 mg/ml)	>240
32	Mesna (50 mg/ml)	>240
33	Mechlorethamine HCL (1mg/ml)	>240
34	Methotrexate (25 mg/ml)	>240
35	Mitomycin C (0.5 mg/ml)	>240
36	Mitoxantrone (2 mg/ml)	>240
37	Oxaliplatin (2 mg/ml)	>240
38	Paclitaxel (6 mg/ml)	>240

	Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
39	Paraplatin (10 mg/ml)	>240
40	Pemetrexed (25 mg/ml)	>240
41	Raltitrexed (0.5 mg/ml)	>240
42	Retrovir (10 mg/ml)	>240
43	Rituximab (10 mg/ml)	>240
44	Temsirolimus (25 mg/ml)	>240
45	ThioTEPA (10 mg/ml)	47.8
46	Topotecan (1 mg/ml)	>240
47	Triclosan (2 mg/ml)	>240
48	Trisenox (1 mg/ml)	>240
49	Vinblastine (1 mg/ml)	>240
50	Vincristine Sulfate (1 mg/ml)	>240
51	Vinorelbine (10 mg/ml)	>240
52	Zoledronic Acid (0.8 mg/ml)	>240

ORDERING INFORMATION - CASES ONLY
Delivery per contract or custom quote upon request

Size	SKU 100 EA/BX 10 BX/CA	SKU 200 EA/BX 10 BX/CA	GTIN - 100 CT GTIN - 200 CT
Small	JVBBCHSMC	JVB200SMC	850020592555 850020592050
Medium	JVBBCHMDC	JVB200MDC	850020592562 850020592180
Large	JVBBCHLGC	JVB200LGC	850020592579 850020592364
X-Large	JVBBCHXLC	JVB200XLC	850020592586 850020592401

For More Information:

855.615.8633

Info@J2MedicalSupply.com

©Medical Supply 2021. All rights reserved.